

RECOVERY IS KEY

Polar Roll-On Massager - GF-POL-PDQ - Display GF-POLROM - Single

- **Stainless Roller Ball** stays ice cold
- **Comfortable, soft rubber grip** fits in your palm
- **Smooth roll-on action**
- **Ice-pack-type, cold rolling massage**
- **Durable TPE and stainless components**
- **Highly portable and compact**
- **Product Size: 3" x 3"**
- **Downloadable Training Poster**



Smooth action,
stainless steel
polar ball swivels
in all directions

Gray retaining ring
unscrews to remove
the ball for cleaning
or freezing



Place in freezer for
approximately
2 hours between uses

