

## PHOTO GEL

Wrist Rests

• Non-skid backing keeps wrist support in place

 Microban® antimicrobial protection keeps product cleaner

 Self adjusting gel support conforms to your wrist for personalized comfort







## **Specifications**

Item #	Description	Color	Covering	Dimensions
9549901	Photo Gel Mouse Pad / Wrist Rest	Black Chevron	Polyurethane	7/8"H x 7 7/8" W x 9 ½" D
9550001	Photo Gel Keyboard Wrist Rest	Black Chevron	Polyurethane	<sup>3</sup> / <sub>4</sub> " H x 18 <sup>9</sup> / <sub>16</sub> " W x 3 <sup>5</sup> / <sub>16</sub> " D

©2016 Fellowes, Inc.

100 YEARS

Since 1917

To receive your check by mail: Clip the "Part Number" from the back cover of the original product manual. Mail with a copy of the receipt, invoice, or packing slip, and the completed form below.

Company*:	Phone*:Email:				
Name & Title*:	Place of Purchase*:				
Address*:	Would you like to participate in market research studies?	☐ Yes	□ No		
City*: State: Zip:	Would you like to receive tips and special offers via e-mail?	☐ Yes	□ No		

Offer valid only on purchases of the Fellowes® Photo Gel Wrist Rests listed in this flyer made between 1/1/17 and 12/31/17. Submit this form along with the original product package UPC and the original invoice dated between 1/1/17 and 12/31/17. Submissions must be postmarked by 1/31/18. Allow 6-8 weeks for delivery. Cannot be combined with any other offer. The check will be issued to the purchaser on the original invoice. Not responsible for lost, late, misdirected, stolen, postage due or undeliverable mail. Illegible or incomplete requests will not be honored. Fellowes dealers, distributors and their employees are ineligible. Void where prohibited, licensed, restricted or taxed by law. Offer only valid for purchases made in the USA. Not valid on clearance items. This offer is limited to consumers; it is not available to clubs or organizations. Offer good in USA only; check by mail paid in US dollars.

Mail to:

Fellowes® Photo Gel Wrist Rests Department 11252 P.O. Box 5017 Stacy, MN 55078-5017